

## Integrating the Healthcare Enterprise



5

# **IHE Patient Care Coordination (PCC) Technical Framework Supplement**

## **Antepartum Profiles (APE, APHP, APL, and APS)**

10

## **Trial Implementation**

15

20 Date: September 9, 2011  
Author: PCC Technical Committee  
Email: [pcc@ihe.net](mailto:pcc@ihe.net)

## Foreword

25 This is a supplement to the IHE Patient Care Coordination Technical Framework V7.0. Each  
supplement undergoes a process of public comment and trial implementation before being  
incorporated into the volumes of the Technical Frameworks.

This supplement is submitted for Trial Implementation as of September 9, 2011 and will be  
available for testing at subsequent IHE Connectathons. The supplement may be amended based  
on the results of testing. Following successful testing it will be incorporated into the PCC  
30 Technical Framework. Comments are invited and can be submitted at  
<http://www.ihe.net/pcc/pcccomments.cfm> or by email to [pcc@ihe.net](mailto:pcc@ihe.net).

This supplement describes changes to the existing technical framework documents and where  
indicated amends text by addition (**bold underline**) or removal (**~~bold strikethrough~~**), as well as  
addition of large new sections introduced by editor’s instructions to “add new text” or similar,  
35 which for readability are not bolded or underlined.

“Boxed” instructions like the sample below indicate to the Volume Editor how to integrate the  
relevant section(s) into the relevant Technical Framework volume:

<i>Replace Section X.X by the following:</i>
--

40

General information about IHE can be found at: [www.ihe.net](http://www.ihe.net)

Information about the IHE Patient Care Coordination domain can be found at:  
<http://www.ihe.net/Domains/index.cfm>

45 Information about the structure of IHE Technical Frameworks and Supplements can be found at:  
<http://www.ihe.net/About/process.cfm> and <http://www.ihe.net/profiles/index.cfm>

The current version of the IHE Technical Framework can be found at:  
[http://www.ihe.net/Technical\\_Framework/index.cfm](http://www.ihe.net/Technical_Framework/index.cfm)

50 **CONTENTS**

	<b>INTRODUCTION.....</b>	<b>5</b>
	HOW TO READ THE ANTEPARTUM PROFILES SUPPLEMENT.....	5
	HOW TO ACCESS THE REFERENCE MATERIAL.....	6
55	OPEN ISSUES AND QUESTIONS.....	6
	CLOSED ISSUES.....	6
	<b>VOLUME 1 – PROFILES .....</b>	<b>7</b>
	1.5 COPYRIGHT PERMISSIONS.....	7
	2.4 DEPENDENCIES OF THE PCC INTEGRATION PROFILES.....	7
60	2.5 HISTORY OF ANNUAL CHANGES .....	7
	<b>W ANTEPARTUM HISTORY AND PHYSICAL CONTENT PROFILE (APHP) .....</b>	<b>8</b>
	W.1 PURPOSE AND SCOPE .....	8
	W.2 PROCESS FLOW .....	8
65	W.2.1 Use Cases.....	8
	W.2.2 Diagrams .....	8
	W.3 ACTORS/TRANSACTIONS.....	8
	W.3.1 Requirements of Actors .....	9
	W.4 OPTIONS.....	9
	W.5 GROUPINGS .....	9
70	W.6 SECURITY CONSIDERATIONS.....	9
	W.7 CONTENT MODULES .....	9
	<b>X ANTEPARTUM SUMMARY CONTENT PROFILE (APS).....</b>	<b>11</b>
	X.1 PURPOSE AND SCOPE.....	11
	X.2 PROCESS FLOW .....	11
75	X.2.1 Use Cases .....	11
	X.2.2 Diagrams.....	11
	X.3 ACTORS/TRANSACTIONS .....	11
	X.3.1 Requirements of Actors.....	12
	X.4 OPTIONS.....	12
80	X.5 GROUPINGS .....	12
	X.6 SECURITY CONSIDERATIONS .....	12
	X.7 CONTENT MODULES.....	12
	<b>Y ANTEPARTUM LABORATORY CONTENT PROFILE (APL).....</b>	<b>14</b>
	Y.1 PURPOSE AND SCOPE.....	14
85	Y.2 PROCESS FLOW.....	14
	Y.2.1 Use Cases .....	14
	Y.2.2 Diagrams.....	14
	Y.3 ACTORS/TRANSACTIONS.....	14
	Y.3.1 Requirements of Actors.....	15
90	Y.4 OPTIONS .....	15
	Y.5 GROUPINGS .....	15

	Y.6 SECURITY CONSIDERATIONS .....	15
	Y.7 CONTENT MODULES.....	15
	<b>Z ANTEPARTUM EDUCATION CONTENT PROFILE (APE) .....</b>	<b>16</b>
95	Z.1 PURPOSE AND SCOPE.....	16
	Z.2 PROCESS FLOW.....	16
	Z.2.1 Use Cases.....	16
	Z.2.2 Diagrams .....	16
100	Z.3 ACTORS/TRANSACTIONS.....	16
	Z.3.1 Requirements of Actors.....	17
	Z.4 OPTIONS .....	17
	Z.5 GROUPINGS .....	17
	Z.6 SECURITY CONSIDERATIONS .....	17
	Z.7 CONTENT MODULES.....	17
105	<b>GLOSSARY .....</b>	<b>19</b>
	<b>VOLUME 2 – TRANSACTIONS AND CONTENT MODULES.....</b>	<b>20</b>
	<b>5 NAMESPACES AND VOCABULARIES.....</b>	<b>21</b>
	5.1 IHE FORMAT CODES.....	21
	<b>6 PCC CONTENT MODULES .....</b>	<b>22</b>
110	6.2 FOLDER CONTENT MODULES.....	22
	6.3 HL7 VERSION 3.0 CONTENT MODULES.....	22
	6.3.1 CDA Document Content Modules.....	22
	6.3.1.A Antepartum History and Physical Specification 1.3.6.1.4.1.19376.1.5.3.1.1.16.1.1 .....	22
115	6.3.1.A.1 Format Code .....	22
	6.3.1.A.2 LOINC Code.....	22
	6.3.1.A.3 Standards.....	22
	6.3.1.A.4 Specification .....	22
	6.3.1.A.5 Conformance .....	24
120	6.3.1.B Antepartum Summary Specification 1.3.6.1.4.1.19376.1.5.3.1.1.11.2 .....	26
	6.3.1.B.1 Format Code .....	26
	6.3.1.B.2 LOINC Code.....	26
	6.3.1.B.3 Standards.....	26
	6.3.1.B.4 Specification .....	27
125	6.3.1.B.5 Conformance .....	29
	6.3.1.C Antepartum Laboratory Specification 1.3.6.1.4.1.19376.1.5.3.1.1.16.1.2 .....	31
	6.3.1.C.1 Format Code .....	31
	6.3.1.C.2 LOINC Code.....	31
	6.3.1.C.3 Standards.....	31
	6.3.1.C.4 Specification.....	31
130	6.3.1.C.5 Conformance .....	31
	6.3.1.D Antepartum Education Specification 1.3.6.1.4.1.19376.1.5.3.1.1.16.1.3 .....	32
	6.3.1.D.1 Format Code.....	32
	6.3.1.D.2 LOINC Code.....	32
	6.3.1.D.3 Standards.....	32
135	6.3.1.D.4 Specification .....	33
	6.3.1.D.5 Conformance .....	33

## Introduction

140 This supplement is written for Trial Implementation. It is written as changes to the documents listed below. The reader should have already read and understood these documents:

1. [PCC Technical Framework Volume 1, Revision 7.0](#)
2. [PCC Technical Framework Volume 2, Revision 7.0](#)

This supplement also references other documents<sup>1</sup>. The reader should have already read and understood these documents:

- 145
1. [IT Infrastructure Technical Framework Volume 1, Revision 8.0](#)
  2. [IT Infrastructure Technical Framework Volume 2, Revision 8.0](#)
  3. [IT Infrastructure Technical Framework Volume 3, Revision 8.0](#)
  4. [The Patient Identifier Cross-Reference \(PIX\) and Patient Demographic Query \(PDQ\) HL7 v3 Supplement to the IT Infrastructure Technical Framework.](#)
- 150
5. HL7 and other standards documents referenced in Volume 1 and Volume 2

## How to read the Antepartum Profiles supplement

This supplement contains 4 content profiles – Antepartum History and Physical (APHP), Antepartum Summary (APS), Antepartum Laboratory (APL), and Antepartum Education (APE). Antepartum Record (APR) is no longer a profile as it has been subsumed by Perinatal Workflow (PW) profile.

Please see the below documents that will need to be referenced to fully understand the profiles in this supplement. Each document has a short description describing what is contained.

1. **Perinatal Workflow (PW):** makes use of the antepartum, labor and delivery, postpartum, and newborn delivery profiles (some are in this supplement and many are in other supplements).
- 155 2. **Content Modules Supplement:** This document contains all PCC Section Templates, Entry Templates and Value Sets that are not in Final Text (that is, they are not in the Technical Framework Volume 2).
3. **PCC Technical Framework Volume 2:** This contains all PCC Section Templates, Entry Templates and Value Sets (among other things) that are in Final Text.

160

---

<sup>1</sup> The first four documents can be located on the IHE Website at [http://www.ihe.net/Technical\\_Framework/index.cfm#IT](http://www.ihe.net/Technical_Framework/index.cfm#IT). The remaining documents can be obtained from their respective publisher.

## How to Access the Reference Material

To access the latest version of the Perinatal Workflow, CDA Content Modules or the PCC Technical Framework, navigate to [http://www.ihe.net/Technical\\_Framework/index.cfm#pcc](http://www.ihe.net/Technical_Framework/index.cfm#pcc)

## Open Issues and Questions

- 165
1. Antepartum Laboratory (APL) currently consists of one document content module – should lab results instead be part of Antepartum Summary as a Coded Results section? This Coded Results section can still point to the Antepartum Laboratory value set as APL currently does now. This approach is consistent with how other data is structured in similar PCC content profiles.

## Closed Issues

170

NA

## Volume 1 – Profiles

*Add the following to section 1.5*

### 175 **1.5 Copyright Permissions**

*Add the following to section 2.4*

### **2.4 Dependencies of the PCC Integration Profiles**

<Profile Name>	<?>	<?>	<->
----------------	-----	-----	-----

180 *Add the following to section 2.5*

### **2.5 History of Annual Changes**

*Add Section W*

## **W Antepartum History and Physical Content Profile (APHP)**

185 Antepartum History and Physical is a content profile that defines the structure of the data that is often collected at the initial ambulatory office visit for a pregnant patient. It includes, but is not limited to demographics, histories, allergies, review of systems, physical examinations, and vital signs.

### **W.1 Purpose and Scope**

190 *Change referenced section numbering when merged into technical framework*

Obstetric patients in labor and admitted to Labor and Delivery must have a complete summary of their antepartum ambulatory care available at the time of admission for optimal care. While this profile represents a portion of that data, and may be included with other antepartum care documents to make up a patient’s Antepartum Record, it may also be created as a standalone document as it is often transmitted separately from those other antepartum care documents. For further information on the full Perinatal Workflow overview and use cases see Perinatal Workflow section X.1.2 and section X.2.

### **W.2 Process Flow**

#### **W.2.1 Use Cases**

*Change referenced section numbering when merged into technical framework*

For applicable use cases see Perinatal Workflow section X.2.1.

#### **W.2.2 Diagrams**

205 *Change referenced section numbering when merged into technical framework*

For applicable diagrams see Perinatal Workflow section X.2.2.

### **W.3 Actors/Transactions**

210 There are two actors in this profile, the Content Creator and the Content Consumer. Content is created by a Content Creator and is to be consumed by a Content Consumer. The sharing or transmission of content from one actor to the other is addressed by the appropriate use of IHE profiles described below, and is out of scope of this profile. A Document Source or a Portable Media Creator may embody the Content Creator Actor. A Document Consumer, a Document Recipient, or a Portable Media Importer may embody the Content Consumer Actor. The sharing or transmission of content or updates from one actor to the other is addressed by the use of



appropriate IHE profiles described in the section on Content Bindings with XDS, XDM and XDR in PCC TF-2:4.1.



220

**Figure X.3-1. Actor Diagram**

### W.3.1 Requirements of Actors

### W.4 Options

Options that may be selected for this Content Profile are listed in the table W.4-1 along with the Actors to which they apply. Dependencies between options when applicable are specified in notes.

**Table W.4-1. Antepartum History and Physical Options**

Actor	Option	Section
Content Consumer	View Option (See Note 1)	PCC TF-2: 3.1.1
	Document Import Option (See Note 1)	PCC TF-2: 3.1.2
	Section Import Option (See Note 1)	PCC TF-2: 3.1.3
	Discrete Data Import Option (See Note 1)	PCC TF-2: 3.1.4
Content Creator	No options defined	

Note 1: The Actor shall support at least one of these options.

230

### W.5 Groupings

### W.6 Security Considerations

### W.7 Content Modules

**Table W.7-1. Antepartum History and Physical Content Modules**

ACOG Antepartum Record Datum	PCC Template Name	PCC Template Id
Demographics	Header Modules	N/A

IHE PCC Technical Framework Supplement – Antepartum History and Physical (APHP),  
Antepartum Summary (APS), Antepartum Laboratory (APL), Antepartum Education (APE)

---

<b>ACOG Antepartum Record Datum</b>	<b>PCC Template Name</b>	<b>PCC Template Id</b>
Chief Complaint	Chief Complaint	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.1
Pregnancy History	Pregnancy History	1.3.6.1.4.1.19376.1.5.3.1.1.5.3.4
Medical History	History of Past Illness	1.3.6.1.4.1.19376.1.5.3.1.3.8
Medical History – Tobacco, Alcohol, Drugs	Coded Social History	1.3.6.1.4.1.19376.1.5.3.1.3.16.1
Medical History – Relevant Family History	Coded Family Medical History	1.3.6.1.4.1.19376.1.5.3.1.3.15
Medications	Medications	1.3.6.1.4.1.19376.1.5.3.1.3.19
Allergies	Allergies and Other Adverse Reactions	1.3.6.1.4.1.19376.1.5.3.1.3.13
Menstrual History/Symptoms Since LMP	Review of Systems	1.3.6.1.4.1.19376.1.5.3.1.3.18
Genetic Screening/Teratology Counseling	Coded Family History Medical History	1.3.6.1.4.1.19376.1.5.3.1.3.15
Infection History	Coded History of Infection	1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1.1
Initial Physical Examination	Coded Physical Exam	1.3.6.1.4.1.19376.1.5.3.1.1.9.15.1
Vital Signs	Coded Vital Signs	1.3.6.1.4.1.19376.1.5.3.1.1.5.3.2
Diagnostic Findings	Coded Results	1.3.6.1.4.1.19376.1.5.3.1.3.28
Surgical History	History of Surgical Procedures	1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2

235

<i>Add Section X</i>
----------------------

## **X Antepartum Summary Content Profile (APS)**

*Change referenced section numbering when merged into technical framework*

240

Antepartum Summary is a content profile that defines the structure for the aggregation of significant events, diagnoses, and plans of care derived from the visits over the course of an antepartum episode. It is represented in part by Estimated Due Dates and a Visit Summary Flowsheet, in which the aggregated data from the ambulatory office visits is recorded, as well as allergies, advance directives, care plans, and selected histories are provided.

245

### **X.1 Purpose and Scope**

*Change referenced section numbering when merged into technical framework*

250

Obstetric patients in labor and admitted to Labor and Delivery must have a complete summary of their antepartum ambulatory care available at the time of admission. While this profile represents a portion of that data and may be included with other antepartum care documents to make up a patient's Antepartum Record, it may also be created as a standalone document and may be transmitted separately from those other antepartum care documents. For further information on the full Perinatal Workflow overview and use cases see Perinatal Workflow section X.1.2 and section X.2.

255

### **X.2 Process Flow**

#### **X.2.1 Use Cases**

*Change referenced section numbering when merged into technical framework*

260

For applicable use cases see Perinatal Workflow section X.2.1.

#### **X.2.2 Diagrams**

*Change referenced section numbering when merged into technical framework*

For applicable diagrams see Perinatal Workflow section X.2.2.

265

### **X.3 Actors/Transactions**

There are two actors in this profile, the Content Creator and the Content Consumer. Content is created by a Content Creator and is to be consumed by a Content Consumer. The sharing or transmission of content from one actor to the other is addressed by the appropriate use of IHE profiles described below, and is out of scope of this profile. A Document Source or a Portable Media Creator may embody the Content Creator Actor. A Document Consumer, a Document

270

Recipient, or a Portable Media Importer may embody the Content Consumer Actor. The sharing or transmission of content or updates from one actor to the other is addressed by the use of appropriate IHE profiles described in the section on Content Bindings with XDS, XDM and XDR in PCC TF-2:4.1.

275



**Figure X.3-1. Actor Diagram**

### X.3.1 Requirements of Actors

280 **X.4 Options**

Options that may be selected for this Content Profile are listed in the table X.4-1 along with the Actors to which they apply. Dependencies between options when applicable are specified in notes.

**Table X.4-1. Antepartum Summary Options**

Actor	Option	Section
Content Consumer	View Option (See Note 1)	PCC TF-2: 3.1.1
	Document Import Option (See Note 1)	PCC TF-2: 3.1.2
	Section Import Option (See Note 1)	PCC TF-2: 3.1.3
	Discrete Data Import Option (See Note 1)	PCC TF-2: 3.1.4
Content Creator	No options defined	

285

Note 1: The Actor shall support at least one of these options.

### X.5 Groupings

### X.6 Security Considerations

290 **X.7 Content Modules**

**Table X.7-1. Antepartum Summary Content Modules**

<b>ACOG Antepartum Record Datum</b>	<b>PCC Template Name</b>	<b>PCC Template Id</b>
Drug Allergy/Latex Allergy	Allergies and Other Adverse Reactions	1.3.6.1.4.1.19376.1.5.3.1.3.13
Is Blood Transfusion Acceptable	Advance Directives	1.3.6.1.4.1.19376.1.5.3.1.3.34
Antepartum Anesthesia Consult Planned	Care Plan	1.3.6.1.4.1.19376.1.5.3.1.3.31
Problems/Plans	Problems	1.3.6.1.4.1.19376.1.5.3.1.3.6
Medication List	Active Medications	1.3.6.1.4.1.19376.1.5.3.1.3.19
EDD Confirmation/18-20 Week EDD Update	Estimated Delivery Dates	1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.1
Prepregnancy Weight	Antepartum Visit Summary Flowsheet	1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.2
Visit Flowsheet	Antepartum Visit Summary Flowsheet	1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.2
	Coded Antenatal Testing and Surveillance	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5.1
	History of Surgical Procedures	1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2

295

*Add Section Y*

## Y Antepartum Laboratory Content Profile (APL)

*Change referenced section numbering when merged into technical framework*

300 Antepartum Laboratory is a content profile that provides the data structure to represent results from standard laboratory tests administered during the antepartum episode.

### Y.1 Purpose and Scope

*Change referenced section numbering when merged into technical framework*

305 Obstetric patients in labor and admitted to Labor and Delivery must have a complete summary of their antepartum ambulatory care available at the time of admission. While this profile represents a portion of that data and may be included with other antepartum care documents to make up a patient's Antepartum Record, it may also be created as a standalone document and may be transmitted separately from those other antepartum care documents. For further information on the full Perinatal Workflow overview and use cases see Perinatal Workflow section X.1.2 and section X.2.

### Y.2 Process Flow

#### Y.2.1 Use Cases

*Change referenced section numbering when merged into technical framework*

315 For applicable use cases see Perinatal Workflow section X.2.1.

#### Y.2.2 Diagrams

*Change referenced section numbering when merged into technical framework*

For applicable diagrams see Perinatal Workflow section X.2.2.

### Y.3 Actors/Transactions

320 There are two actors in this profile, the Content Creator and the Content Consumer. Content is created by a Content Creator and is to be consumed by a Content Consumer. The sharing or transmission of content from one actor to the other is addressed by the appropriate use of IHE profiles described below, and is out of scope of this profile. A Document Source or a Portable Media Creator may embody the Content Creator Actor. A Document Consumer, a Document Recipient, or a Portable Media Importer may embody the Content Consumer Actor. The sharing or transmission of content or updates from one actor to the other is addressed by the use of appropriate IHE profiles described in the section on Content Bindings with XDS, XDM and XDR in PCC TF-2:4.1.



330

**Figure Y.3-1. Actor Diagram**

### Y.3.1 Requirements of Actors

### Y.4 Options

335 Options that may be selected for this Content Profile are listed in the table Y.4-1 along with the Actors to which they apply. Dependencies between options when applicable are specified in notes.

**Table Y.4-1. Antepartum Laboratory Options**

Actor	Option	Section
Content Consumer	View Option (See Note 1) Document Import Option (See Note 1) Section Import Option (See Note 1) Discrete Data Import Option (See Note 1)	PCC TF-2: 3.1.1 PCC TF-2: 3.1.2 PCC TF-2: 3.1.3 PCC TF-2: 3.1.4
Content Creator	No options defined	

Note 1: The Actor shall support at least one of these options.

340

### Y.5 Groupings

### Y.6 Security Considerations

### Y.7 Content Modules

345 There are no content module mappings for Antepartum Laboratory as XD-LAB profile is referenced in full.

*Add Section Z*

## **Z Antepartum Education Content Profile (APE)**

*Change referenced section numbering when merged into technical framework*

350

Antepartum Education is a content profile that provides the document structure to represent educational material provided during the office visit(s) for the antepartum episode.

### **Z.1 Purpose and Scope**

*Change referenced section numbering when merged into technical framework*

355

Obstetric patients in labor and admitted to Labor and Delivery must have a complete summary of their antepartum ambulatory care available at the time of admission. While this profile represents a portion of that data and may be included with other antepartum care documents to make up a patient's Antepartum Record, it may also be created as a standalone document and may be transmitted separately from those other antepartum care documents. For further information on the full Perinatal Workflow overview and use cases see Perinatal Workflow section X.1.2 and section X.2.

360

### **Z.2 Process Flow**

#### **Z.2.1 Use Cases**

365

*Change referenced section numbering when merged into technical framework*

For applicable use cases see Perinatal Workflow section X.2.1.

#### **Z.2.2 Diagrams**

*Change referenced section numbering when merged into technical framework*

370

For applicable diagrams see Perinatal Workflow section X.2.2.

### **Z.3 Actors/Transactions**

There are two actors in this profile, the Content Creator and the Content Consumer. Content is created by a Content Creator and is to be consumed by a Content Consumer. The sharing or transmission of content from one actor to the other is addressed by the appropriate use of IHE profiles described below, and is out of scope of this profile. A Document Source or a Portable Media Creator may embody the Content Creator Actor. A Document Consumer, a Document Recipient, or a Portable Media Importer may embody the Content Consumer Actor. The sharing or transmission of content or updates from one actor to the other is addressed by the use of

375



380 appropriate IHE profiles described in the section on Content Bindings with XDS, XDM and  
 XDR in PCC TF-2:4.1.



Figure Z.3-1. Actor Diagram

385

### Z.3.1 Requirements of Actors

### Z.4 Options

Options that may be selected for this Content Profile are listed in the table Z.4-1 along with the  
 Actors to which they apply. Dependencies between options when applicable are specified in  
 notes.

390

Table Z.4-1. Antepartum Education Options

Actor	Option	Section
Content Consumer	View Option (See Note 1)	PCC TF-2: 3.1.1
	Document Import Option (See Note 1)	PCC TF-2: 3.1.2
	Section Import Option (See Note 1)	PCC TF-2: 3.1.3
	Discrete Data Import Option (See Note 1)	PCC TF-2: 3.1.4
Content Creator	No options defined	

Note 1: The Actor shall support at least one of these options.

### Z.5 Groupings

395 **Z.6 Security Considerations**

### Z.7 Content Modules

**Table Z.7-1. Antepartum Education Content Modules**

<b>ACOG Antepartum Record Datum</b>	<b>PCC Template Name</b>	<b>PCC Template Id</b>
Education, First Trimester	Patient Education	1.3.6.1.4.1.19376.1.5.3.1.1.9.38
Education, Second Trimester	Patient Education	1.3.6.1.4.1.19376.1.5.3.1.1.9.38
Education, Third Trimester	Patient Education	1.3.6.1.4.1.19376.1.5.3.1.1.9.38

400

## **Glossary**

*Add the following terms to the Glossary:*

405

## **Volume 2 – Transactions and Content Modules**

## 5 Namespaces and Vocabularies

codeSystem	codeSystemName	Description

410

### 5.1 IHE Format Codes

Profile	Format Code	Media Type	Template ID
Antepartum History and Physical (APHP)	urn:ihe:pcc:aphp:2008	text/xml	1.3.6.1.4.1.19376.1.5.3.1.1.16.1.1
Antepartum Summary (APS)	urn:ihe:pcc:aps:2007	text/xml	1.3.6.1.4.1.19376.1.5.3.1.1.11.2
Antepartum Laboratory (APL)	urn:ihe:pcc:apl:2008	text/xml	1.3.6.1.4.1.19376.1.5.3.1.1.16.1.2
Antepartum Education (APE)	urn:ihe:pcc:ape:2008	text/xml	1.3.6.1.4.1.19376.1.5.3.1.1.16.1.3

## 6 PCC Content Modules

### 415 6.2 Folder Content Modules

See Perinatal Workflow section 6.2.A.

### 6.3 HL7 Version 3.0 Content Modules

#### 6.3.1 CDA Document Content Modules

*Add section 6.3.1.A*

#### 420 6.3.1.A Antepartum History and Physical Specification 1.3.6.1.4.1.19376.1.5.3.1.1.16.1.1

The Antepartum History and Physical contains a record of the History and Physical usually performed during the initial visit. This document content module is a History and Physical and inherits all header constraints from History and Physical (1.3.6.1.4.1.19376.1.5.3.1.1.16.1.4)

#### 425 6.3.1.A.1 Format Code

The XDSDocumentEntry format code for this content is **urn:ihe:pcc:aphp:2008**

#### 6.3.1.A.2 LOINC Code

The LOINC code for this document is **34117-2 HISTORY AND PHYSICAL**

#### 6.3.1.A.3 Standards

<b>CCD</b>	ASTM/HL7 Continuity of Care Document
<b>CDAR2</b>	<a href="#">HL7 CDA Release 2.0</a>
<b>ACOG AR</b>	<a href="#">American College of Obstetricians and Gynecologists (ACOG), Antepartum Record</a>
<b>LOINC</b>	<a href="#">Logical Observation Identifiers, Names and Codes</a>
<b>SNOMED</b>	<a href="#">Systemized Nomenclature for Medicine</a>
<b>CDTHP</b>	<a href="#">CDA for Common Document Types History and Physical Notes (DSTU)</a>

#### 430 6.3.1.A.4 Specification

This section references content modules using Template Id as the key identifier. Definitions of the modules are found in either:

- IHE Patient Care Coordination Volume 2: Final Text
- IHE PCC Content Modules 2011 Supplement

435

**Table 6.3.1.A.4-1. Antepartum History and Physical Specification**

Template Name	Opt	Section Template Id / Location	Value Set Template Id
Spouse	R	1.3.6.1.4.1.19376.1.5.3.1.2.4.1 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.2.5	N/A
Natural Father of Fetus	R	1.3.6.1.4.1.19376.1.5.3.1.2.4.2 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.2.6	N/A
Ethnicity The ethnicity of the patient should be recorded	R2		N/A
Chief Complaint	R	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.1 IHE PCC 2:6.3.3.1.3	N/A
History of Present Illness	R	1.3.6.1.4.1.19376.1.5.3.1.3.4 IHE PCC 2:6.3.3.1.3	N/A
History of Past Illness This section is the same as it is for History and Physical, however it SHALL contain entries and SHOULD use codes as specified in the Antepartum History and Physical History of Past Illness Value Set. A negative diagnosis SHALL be recorded with the use of the negation indicator attribute. If the data is not present or not available within the system no entry is required.	R	1.3.6.1.4.1.19376.1.5.3.1.3.8 IHE PCC 2:6.3.3.2.1	1.3.6.1.4.1.19376.1.5.3.1.1.16.5.1
Coded History of Infection This section SHOULD use the codes specified in the Antepartum History of Infection Value Set.	R	1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1.1 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.3.2.37	1.3.6.1.4.1.19376.1.5.3.1.1.16.5.6
Pregnancy History	R	1.3.6.1.4.1.19376.1.5.3.1.1.5.3.4 IHE PCC 2:6.3.3.2.18	N/A
Coded Social History	R	1.3.6.1.4.1.19376.1.5.3.1.3.16.1 PCC TF Supplement CDA Content Modules (TI) Vol 2:6.3.3.2.36	N/A
Coded Family Medical History This section is the same as it is for History & Physical, however it SHALL contain Genetic Screening and Teratology Counseling information as specified in the	R	1.3.6.1.4.1.19376.1.5.3.1.3.15 IHE PCC 2:6.3.3.2.13	1.3.6.1.4.1.19376.1.5.3.1.1.16.5.4

Template Name	Opt	Section Template Id / Location	Value Set Template Id
Antepartum Family History and Genetic Screening Value Set. If the data is not present or not available within the system no entries are required.			
Allergies and Other Adverse Reactions	R	1.3.6.1.4.1.19376.1.5.3.1.3.13 IHE PCC 2:6.3.3.2.11	N/A
Review of Systems This section is the same as it is for History & Physical, however it SHALL include organizers for Menstrual History and MAY include entries for general review of systems data. The Menstrual History entries SHOULD use the codes specified in the Antepartum Menstrual History Value Set. The section code value for the Menstrual History organizer SHALL be '49033-4'. The codeSystem is '2.16.840.1.113883.6.1'. The codeSystemName is 'LOINC'.	R	1.3.6.1.4.1.19376.1.5.3.1.3.18 IHE PCC 2:6.3.3.2.16	1.3.6.1.4.1.19376.1.5.3.1.1.16.5.5
Coded Physical Exam This section SHALL include Coded Vital Signs (1.3.6.1.4.1.19376.1.5.3.1.1.5.3.2) as a subsection.	R	1.3.6.1.4.1.19376.1.5.3.1.1.9.15.1 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.3.4.30	N/A
History of Surgical Procedures	R2	1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.3.2.44	N/A

### 6.3.1.A.5 Conformance

440 CDA Release 2.0 documents that conform to the requirements of this content module shall indicate their conformance by the inclusion of the appropriate <templateId> elements in the header of the document. This is shown in the sample document below. A CDA Document may conform to more than one template. This content module inherits from the History and Physical content module, and so must conform to the requirements of that template as well, thus all <templateId> elements shown in the example below shall be included.



```
445 <ClinicalDocument xmlns='urn:hl7-org:v3'>
  <typeId extension="POCD_HD000040" root="2.16.840.1.113883.1.3"/>
  <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.2'/><!--Medical Summary-->
  <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.1.4'/><!--History and Physical-->
  <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.1.1'/><!--Antepartum History and Physical-->
  <id root=' ' extension=' '/>
450 <code code='34117-2' displayName='HISTORY AND PHYSICAL'
  codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'/>
  <title>Antepartum History and Physical</title>
  <effectiveTime value='20080601012005'/>
  <confidentialityCode code='N' displayName='Normal'
455   codeSystem='2.16.840.1.113883.5.25' codeSystemName='Confidentiality' />
  <languageCode code='en-US'/>
  :
  <component><structuredBody>
460   <component>
     <section>
       <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.1'/>
       <!-- Required Chief Complaint Section content -->
     </section>
   </component>
465   <component>
     <section>
       <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.4'/>
       <!-- Required History of Present Illness Section content -->
     </section>
   </component>
470   <component>
     <section>
       <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.8'/>
       <!-- Required History of Past Illness Section content -->
475     </section>
   </component>
   <component>
     <section>
       <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1.1'/>
       <!-- Required Coded History of Infection Section content -->
480     </section>
   </component>
   <component>
     <section>
       <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.5.3.4'/>
       <!-- Required Pregnancy History Section content -->
485     </section>
   </component>
   <component>
     <section>
       <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.16'/>
       <!-- Required Coded Social History Section content -->
490     </section>
   </component>
   <component>
     <section>
       <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.15'/>
       <!-- Required Coded Family Medical History Section content -->
495     </section>
   </component>
   <component>
     <section>
       <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.13'/>
       <!-- Required Allergies and Other Adverse Reactions Section content -->
500     </section>
   </component>
   <component>
     <section>
505     </section>
   </component>
```

```

510     <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.18' />
        <!-- Required Review of Systems Section content -->
    </section>
</component>
<component>
    <section>
515     <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.15' />
        <!-- Required Physical Examination Section content -->
    </section>
</component>
<component>
520     <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.25' />
        <!-- Conditional Vital Signs Section content -->
    </section>
</component>
<component>
525     <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2' />
        <!-- Required if known History of Surgical Procedures Section content -->
    </section>
</component>
530 </structuredBody></component>
</ClinicalDocument>

```

**Figure 6.3.1.A-5-1 Sample Antepartum History and Physical Document**

*Add section 6.3.1.B*

**535 6.3.1.B Antepartum Summary Specification 1.3.6.1.4.1.19376.1.5.3.1.1.11.2**

The Antepartum Summary represents a summary of the most critical information to an antepartum care provider regarding the status of a patient’s pregnancy. This document content module is a Medical Summary and inherits all header constraints from Medical Summary (1.3.6.1.4.1.19376.1.5.3.1.1.2).

**540 6.3.1.B.1 Format Code**

The XDSDocumentEntry format code for this content is **urn:ihe:pcc:aps:2007**

**6.3.1.B.2 LOINC Code**

The LOINC code for this document is **57055-6** Antepartum summary

**6.3.1.B.3 Standards**

<b>CCD</b>	ASTM/HL7 Continuity of Care Document
<b>CDAR2</b>	<a href="#">HL7 CDA Release 2.0</a>
<b>ACOG AR</b>	<a href="#">American College of Obstetricians and Gynecologists (ACOG), Antepartum Record</a>
<b>LOINC</b>	<a href="#">Logical Observation Identifiers, Names and Codes</a>
<b>SNOMED</b>	<a href="#">Systemized Nomenclature for Medicine</a>
<b>CDTHP</b>	<a href="#">CDA for Common Document Types History and Physical Notes (DSTU)</a>

545 **6.3.1.B.4 Specification**

This section references content modules using Template ID as the key identifier. Definitions of the modules are found in either:

- IHE Patient Care Coordination Volume 2: Final Text
- IHE PCC Content Modules 2011 Supplement

550

**Table 6.3.1.B.4-1. Antepartum Summary Specification**

Template Name	Opt	Section Template Id / Location	Value Set Template Id
<p>Allergies and Other Adverse Reactions</p> <p>This section is the same as for Medical Summary, however it SHALL include one observation of Latex Allergy which may be negated through the negationInd attribute. Latex Allergy is particularly relevant for Obstetrics because of the frequency of vaginal exams that might involve the use of latex gloves. The observation value code for Latex Allergy is '300916003'. The codeSystem is '2.16.840.1.113883.6.96'. The codeSystemName is 'SNOMED CT'</p>	R	<p>1.3.6.1.4.1.19376.1.5.3.1.3.13</p> <p>IHE PCC 2:6.3.3.2.11</p>	
<p>Advance Directives</p> <p>APS includes an explicit check of patients preference for blood transfusion because the risk of massive hemorrhage during delivery is much higher. This observation SHALL be recorded in the Advance Directives section. A simple observation of "blood transfusion acceptable?" SHALL be included with the value '(xx-bld-transf-ok)'. The codeSystem is '2.16.840.1.113883.6.1'. The codeSystemName is 'LOINC'</p>	R	<p>1.3.6.1.4.1.19376.1.5.3.1.3.34</p> <p>IHE PCC 2:6.3.3.6.5</p>	
<p>Care Plan</p> <p>APS forms SHOULD include an observation stating if an anesthesia consult is planned. When present, the observation value for this observation is '(xx-anest-cons-pland)'. The codeSystem is '2.16.840.1.113883.6.1'. The codeSystemName is 'LOINC'. If the type of anesthesia planned is known, systems SHOULD include</p>	R	<p>1.3.6.1.4.1.19376.1.5.3.1.3.31</p> <p>IHE PCC 2:6.3.3.6.1</p>	

IHE PCC Technical Framework Supplement – Antepartum History and Physical (APHP),  
Antepartum Summary (APS), Antepartum Laboratory (APL), Antepartum Education (APE)

Template Name	Opt	Section Template Id / Location	Value Set Template Id
an observation to represent that data using the LOINC code '(xx-type-of-anesth-pland)' with a CD value including one of the following values: ( General   Epidural   Spinal ) or a Null flavor to represent unknown or not listed.			
Medications Medications should include start and stop date if known.	R	1.3.6.1.4.1.19376.1.5.3.1.3.19 IHE PCC 2:6.3.3.1	
Active Problems Related Plans should be included in the Care Plan section.	R	1.3.6.1.4.1.19376.1.5.3.1.3.6 IHE PCC 2:6.3.3.2.3	
Estimated Delivery Date	R	1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.1 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.3.2.28	
Antepartum Visit Summary Flowsheet	R	1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.2 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.3.9.3	
History of Surgical Procedures	R2	1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.3.2.44	
Coded Antenatal Testing and Surveillance	R2	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5.1 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.3.5.7	

Note: The Antepartum Summary is typically used as a 'living document' where the latest information is added to the end of the flowsheet at each visit. This is different than a typical Medical Summary which typically would not share information until document is complete. Although this pattern of updates is not prohibited by Medical Summary, it is also not typical. For APS documents may be published at the end of each visit, but subsequent updates with a pregnancy SHALL be represented as document replacement by including a <relatedDocument typeCode='REPL'> element as below.

555

560

565

```
<ClinicalDocument xmlns='urn:hl7-org:v3'>
:
<relatedDocument typeCode='REPL'>
<parentDocument>
<id root=' ' extension=' ' />
</parentDocument>
</relatedDocument>
:
</ClinicalDocument>
```

Figure 6.3.1.B.4-1

### **6.3.1.B.5 Conformance**

570 CDA Release 2.0 documents that conform to the requirements of this content module shall  
indicate their conformance by the inclusion of the appropriate <templateId> elements in the  
header of the document. This is shown in the sample document below. A CDA Document may  
conform to more than one template. This content module inherits from the Medical Summary  
content module, and so must conform to the requirements of that template as well, thus all  
575 <templateId> elements shown in the example below shall be included.

```
580 <ClinicalDocument xmlns='urn:hl7-org:v3'>
  <typeId extension="POCD_HD000040" root="2.16.840.1.113883.1.3"/>
  <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.2'/><!--Medical Summary-->
  <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2'/><!--Antepartum Summary-->
  <id root=' ' extension=' '/>
  <code code='57055-6' displayName='Antepartum summary'
    codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
  <title>Antepartum Summary</title>
  <effectiveTime value='20080601012005'/>
  <confidentialityCode code='N' displayName='Normal'
    codeSystem='2.16.840.1.113883.5.25' codeSystemName='Confidentiality' />
  <languageCode code='en-US'/>
  :
  <component><structuredBody>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.13'/>
        <!-- Required Allergies and Other Adverse Reactions Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.34'/>
        <!-- Required Advance Directives Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.31'/>
        <!-- Required Care Plan Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.19'/>
        <!-- Required Medications Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.6'/>
        <!-- Required Problems Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.1'/>
        <!-- Required Estimated Delivery Dates Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.2'/>
        <!-- Required Antepartum Visit Summary Flowsheet Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2'/>
        <!-- Required History of Surgical Procedures Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5.1'/>

```

645

```

        <!-- Required Coded Antenatal Testing and Surveillance Section content -->
        </section>
    </component>
</structuredBody></component>
</ClinicalDocument>
    
```

**Figure 6.3.1.B.5-1. Sample Antepartum Summary Document**

*Add Section 6.3.1.C: Antepartum Lab Specification*

**6.3.1.C Antepartum Laboratory Specification 1.3.6.1.4.1.19376.1.5.3.1.1.16.1.2**

650 The Antepartum Laboratory SHALL follow all constraints as defined in the XD-LAB profile, as described in [LAB TF-3:4](#). There is a suggested code list provided in the Antepartum Laboratory Value Set. Due to the variation possible in these laboratory results, and the potential for new codes representing new types of laboratory data, a tightly constrained code list is not provided.

**6.3.1.C.1 Format Code**

655 The XDSDocumentEntry format code for this content is **urn:ihe:pcc:apl:2008**

**6.3.1.C.2 LOINC Code**

The LOINC code for this document is **26436-6** Laboratory Studies

**6.3.1.C.3 Standards**

<b>CCD</b>	ASTM/HL7 Continuity of Care Document
<b>CDAR2</b>	<a href="#">HL7 CDA Release 2.0</a>
<b>ACOG AR</b>	<a href="#">American College of Obstetricians and Gynecologists (ACOG), Antepartum Record</a>
<b>LOINC</b>	<a href="#">Logical Observation Identifiers, Names and Codes</a>
<b>CDTHP</b>	<a href="#">CDA for Common Document Types History and Physical Notes (DSTU)</a>

**6.3.1.C.4 Specification**

660

**Table 6.3.1.C.4-1. Antepartum Laboratory Specification**

Template Name	Opt	Section Template Id	Value Set Template Id
Laboratory Results The Antepartum Laboratory Value Set should be used to represent laboratory results.	R	N/A	1.3.6.1.4.1.19376.1.5.3.1.1.16.5.7

**6.3.1.C.5 Conformance**

CDA Release 2.0 documents that conform to the requirements of this content module shall indicate their conformance by the inclusion of the appropriate <templateId> elements in the

header of the document. This is shown in the sample document below. A CDA Document may  
665 conform to more than one template. This content module inherits from the XD Lab Report  
content module, and so must conform to the requirements of that template as well, thus all  
<templateId> elements shown in the example below shall be included.

```

670 <ClinicalDocument xmlns='urn:hl7-org:v3'>
    <typeId extension="POCD_HD000040" root="2.16.840.1.113883.1.3"/>
    <templateId root='1.3.6.1.4.1.19376.1.3.3'/>
    <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.1.2'/><!--Antepartum Laboratory-->
    <id root=' ' extension=' '/>
675 <code code='26436-6' displayName='Laboratory Studies'
    codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'/>
    <title>Antepartum Laboratory</title>
    <effectiveTime value='20080601012005'/>
    <confidentialityCode code='N' displayName='Normal'
680 codeSystem='2.16.840.1.113883.5.25' codeSystemName='Confidentiality' />
    <languageCode code='en-US'/>
    :
    <component><structuredBody>
685 </structuredBody></component>
</ClinicalDocument>

```

**Figure 6.3.1.C.5-1 Sample Antepartum Laboratory Document**

*Add Section 6.3.1.D: Antepartum Education Specification*

**6.3.1.D Antepartum Education Specification 1.3.6.1.4.1.19376.1.5.3.1.1.16.1.3**

690 The Antepartum Education contains a list of patient education activities that have occurred or  
have been planned to review with the patient.

**6.3.1.D.1 Format Code**

The XDSDocumentEntry format code for this content is **urn:ihe:pcc:ape:2008**

**6.3.1.D.2 LOINC Code**

695 The LOINC code for this document is **34895-3** EDUCATION NOTE

**6.3.1.D.3 Standards**

<b>CCD</b>	ASTM/HL7 Continuity of Care Document
<b>CDAR2</b>	<a href="#">HL7 CDA Release 2.0</a>
<b>ACOG AR</b>	<a href="#">American College of Obstetricians and Gynecologists (ACOG), Antepartum Record</a>
<b>LOINC</b>	<a href="#">Logical Observation Identifiers, Names and Codes</a>
<b>SNOMED</b>	<a href="#">Systemized Nomenclature for Medicine</a>
<b>CDTHP</b>	<a href="#">CDA for Common Document Types History and Physical Notes (DSTU)</a>



### 6.3.1.D.4 Specification

**Table 6.3.1.D.4-1. Antepartum Education Specification**

Template Name	Opt	Section Template Id	Value Set Template Id
Patient Education This section SHALL follow all constraints as listed in the Patient Education Section, and SHOULD use the codes specified in the Antepartum Education Value Set.	R	1.3.6.1.4.1.19376.1.5.3.1.1.9.38 PCC TF Supplement CDA Content Modules (TI) Vol 2: 6.3.3.6.14	1.3.6.1.4.1.19376.1.5.3.1.1.16.5.8

### 6.3.1.D.5 Conformance

700 CDA Release 2.0 documents that conform to the requirements of this content module shall  
 indicate their conformance by the inclusion of the appropriate <templateId> elements in the  
 header of the document. This is shown in the sample document below. A CDA Document may  
 conform to more than one template. This content module inherits from the Medical Documents  
 content module, and so must conform to the requirements of that template as well, thus all  
 705 <templateId> elements shown in the example below shall be included.

```

710 <ClinicalDocument xmlns='urn:hl7-org:v3'>
    <typeId extension="POCD_HD000040" root="2.16.840.1.113883.1.3"/>
    <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.1'/><!--Medical Document-->
    <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.1.3'/><!--Antepartum Education-->
    <id root=' ' extension=' '/>
    <code code='34895-3' displayName='EDUCATION NOTE'
    codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'/>
715 <title>Antepartum Education</title>
    <effectiveTime value='20080601012005'/>
    <confidentialityCode code='N' displayName='Normal'
    codeSystem='2.16.840.1.113883.5.25' codeSystemName='Confidentiality' />
    <languageCode code='en-US'/>
    :
720 <component><structuredBody>
    <component>
    <section>
    <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.38'/>
    <!-- Required Patient Education Section content -->
725 </section>
    </component>
    </structuredBody></component>
    </ClinicalDocument>
    
```

**Figure 6.3.1.D.5-1. Sample Antepartum Education Document**

730